

Maxson®




# SHARPS CONTAINERS

# Sharps Containers

Occupational Safety is important for us, and clinical practitioners should not be exposed to unnecessary risks. Ranging from small needles to surgical sharps instruments collections, all of our Sharps Containers products conform to major international safety standards minimizing the chance of accidental sharps injury. Each has a temporary and a permanent locking mechanism, and made of puncture resistant materials. Each product has a label with clear fill line and installation procedure on it.



## 3 in 1 Locking System



## Features:



### Dual locking system -

Permanent and temporary, to prevent loss of content during transit or usage.



### Needle detaching holes-

3 needle detaching holes for different syringes: lure lock, slip tip and eccentric tip. Prevent direct contact with needle. Protect the healthcare practitioner.







### International Compliance -

BS 7320: 1900  
OSHA 29 CFR 1910.1030  
ASTM F2132 - 01(2008)e1  
AS 4031: 1992










## Sharps Containers:

	Pocket Size		Small Size		
					
Item No.	MI279	MI207	MI155	MI312	MI235
Volume (L)	0.07	0.6	1	1	1.4
Outer Dimensions (mm)	39 x 39 x 166	100 x 55 x 200	110 x 95 x 190	118 x 115 x 176	110 x 95 x 270
Opening Dimensions (mm)	ø28	37 x 22	60 x 35	260 x 210	60 x 35
Temporary Closure	1	2	-	-	-

	Medium Size			Large Size
				
Item No.	MI128	MI323	MI313	MI218
Volume (L)	5	6	7	40
Outer Dimensions (mm)	210 x 155 x 255	200 x 200 x 290	245 x 175 x 285	352 x 352 x 600
Opening Dimensions (mm)	83 x 33	240 x 380	350 x 122	100 x 300
Temporary Closure	1	-	-	1

## Accessories:

0.6L Paper Stand	1/1.4L Stand	1L Hanger	5L Hanger	7L Hanger	40L Wall Mount	40L Trolley
						

**Maxson®**

### ALLEVA MEDICAL LTD.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong  
Tel: +852 3166 7200 Fax: +852 2355 7663

### ALLEVA MEDICAL USA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403  
Tel: 1 267 428 2487 info@allevamedical.com



[www.allevamedical.com](http://www.allevamedical.com)

ALL products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd.  
Alleva Medical USA LLC, Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks manufactured by Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, Europe, China and other countries. Copyright 2021.

All dimensions, weights, colours and contents may vary.