



CDD100 Dry Seal Chest Drainage Unit Distributor Sell Sheet

CDD100 Dry Seal Chest Drainage Unit





Vide®'s Dry Seal is a 2000mL sterile disposable three-chamber system that offers the reliability of a traditional water seal drain and the convenience of a water-less operation. With a proprietary dry suction control technology, suction can be adjusted easily and precisely with one hand. Featuring needleless sampling port and relief valves, the dry seal drain is built entirely with patient safety and ergonomics in mind.

Indications for use: pneumothorax, hemothorax, mediastinal shift and malignant effusions.


Features:

- ✓ **Colored Liquids For Water Seal**
Pre-filled syringe stored safely at the back of the drain and is ready for use



- ✓ **Suction and Negative Pressure Indicators**
Allow easy visualization of the status of the unit




✓ **Intuitive Sliding Suction Control panel** - no need for finding and filling sterile water!


- ✓ **Negative Pressure Relief Valve**
release automatically during prolonged episode of negative intrapleural pressure

✓ **7-Scale Air Leak Monitor**
for more precise air leak quantification


- ✓ **Vide®'s Calibrated** Dry Suction Chest Drain can maintain consistent suction control at varying negative pressure or flow rate from the source

✓ **Y-Shape Sampling Port**
with self-sampling needless feature, allowing access to fresh blood samples


- ✓ **Knock-over nozzle** to ensure function is retained even after unit is accidentally knocked over

✓ **Low Positive Pressure Relief Valve**
release automatically at 2cmH₂O
- ✓ **8-Discrete Suction Levels** with Suction Control Indicator
5cmH₂O to 40cmH₂O

Accessories:

HC063-0001

Chest Drain Basket
with Clamp & Steel Hanger



HC063-0002

Chest Drain Basket with Clamp



ST380-1031

Small-bore Adaptor (22-26 Fr)



Shipping / Packaging information:

| Item No. | Type | Volume(mL) | Packaging |
|----------|------------------------------|------------|----------------|
| CDD100 | Dry Seal Chest Drainage Unit | 2000mL | 5pcs / Carton; |



ALLEVA MEDICAL LTD.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong
Tel: +852 3166 7200 Fax: +852 2355 7663

ALLEVA MEDICAL USA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403
Tel: 1 267 428 2487 info@allevamedical.com



www.allevamedical.com

ALL products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd.
Alleva Medical USA LLC, Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks manufactured by Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, Europe, China and other countries. Copyright 2021.

All dimensions, weights, colours and contents may vary.