



ST380 Three-Chamber Water Seal Chest Drain Distributor Sell Sheet

ST380 Three-Chamber Water Seal Chest Drain



Vide's Water Seal Chest Drainage Unit is a 2,000ml sterile disposable three-chamber chest drain. It is designed as a single collection unit used to remove air and fluid from a patient's intrapleural cavity, as a result of pneumothorax, hemothorax, mediastinal shift or other complications that affect the normal breathing resulting from a loss in intrapleural negative pressure.

Safe • Sturdy • Ergonomic

The ST380 chest drain can be used as a gravity drain or with a vacuum source.

Features:



Suction Control (A):

Allows for a maximum 20cmH₂O suction.
Water level adjustable through self-sealing port.



Diagnostics with Air Leak Monitor (Water Seal B):

Allows for a maximum 20cmH₂O suction.
Water level adjustable through self-sealing port.



Self Sealing Sampling Port (C):

Direct and fresh sample extraction by 24Ga
needle/syringe through self-sealing port.



Reliable 3-Column Collection Chambers:

- 1ml increments up to 100ml
- 2ml increments up to 220ml
- 10ml increments up to 2000ml
- Patient Name and Time can be recorded directly on the unit



SAFE

Easy to Read Front Panel
ISO 10079-3 Compliant



IMPACT RESISTANT

Sturdy Construction With
Hanger and Basket
Size: 24.5*36.0*6.0cm



PATIENT PROTECTION

Automatic Release Valves for both
Positive & Negative Pressure.
Manual Release
Valve Override.

Safety Features:



- ✓ **Low Positive Pressure Relief Valve** opens automatically at 4cm H₂O



- ✓ **High Negativity Relief Valve** permits manual venting of high negative pressure, offering maximum patient safety.



- ✓ Kink Free Tubing Design
- ✓ Tube Management
- ✓ Bed Side Hanger



- ✓ Graduated Cup
- ✓ Control Valve
- ✓ Floor Stand

Accessories:

HC063-0001

Chest Drain Basket
with Clamp & Steel Hanger



HC063-0002

Chest Drain Basket with Clamp



ST380-1031

Small-bore Adaptor (22-36 Fr)



Shipping / Packaging information:

Order code	Product description	Volume(mL)	Packaging
ST380	3-chamber water seal and suction control	2000mL	5pcs / Carton;

ALL products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd.
Alleva Medical USA LLC, Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks manufactured by Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, Europe, China and other countries. Copyright 2021.

All dimensions, weights, colours and contents may vary.

 **ALLEVA** MEDICAL LTD.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong
Tel: +852 3166 7200 Fax: +852 2355 7663

 **ALLEVA** MEDICAL USA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403
Tel: 1 267 428 2487 info@allevamedical.com



www.allevamedical.com